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A second upper jaw is now visible in the list on the left.

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Now the upper jaw can be scanned in the usual way.

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After completing the scan, the target region (bite fork carrier) must be scanned, for this purpose the next upper jaw scan is selected by clicking on "OK".

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! It is important to record from the occlusion to the landing pad on the sagittal rod !

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Now place the transfer sheet with the AxioPrisa bite fork carrier on the patient.

Instructions for using the AxioQuick transfer sheet can be found here: https://www.sam-dental.de/downloads/

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Put a small amount of silicone into the recess of the AxioJig.

Führe nun den AxioJig an die Frontzähne. Now guide the AxioJig to the anterior teeth.

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The patient must now bite into the silicone and fix the AxioJig with the lower jaw. When the silicone is hard, close the screw on the right side of the sagittal bar and the patient can open the mouth. Now the AxioPrisa can be tensioned.

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Next you need to complete the mandible ...

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Wait until the optimisation is complete.

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The "Export" button on the right-hand side can be used to save the files on the local PC so that they can then be aligned in the AxioPrisa software.

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! Important note on processing scans in AxioSnapMount !

Due to possible artifacts or errors in the scan, the generation process may show an error. To prevent this, it is important that the scan of the maxilla as well as the mandible is therefore freed from errors (cleaning the scans in the scan software / CAD software).

Holes in the scan or superimposed surfaces are susceptible.

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3disc.com/contact

For questions about 3DISC hardware or software please contact 3DISC Support:

All enquiries can be made conveniently via the online Help Centre:

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Info@smart-dental.de

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it-support@sam-dental.de

For questions about Axioprisa hardware or software contact us at SAM Präzisionstechnik:

For general questions about the product:

Info@smart-dental.de

In case of technical problems:

it-support@sam-dental.de